

ChIP-Chip

by Kengo Morohashi (morohashi.1@osu.edu)

1. Materials

1.1 Chromatin-IP

See protocol Chromatin-IP

1.2 Amplification of Immunoprecipitated DNA Targets

5X sequenase buffer (USB#70775)

Sequenase (13 U/ul) (USB #70775Y)

primer A (GTTTCCCAGTCACGATCNNNNNNNNN, HPLC purified)

primer B (GTTTCCCAGTCACGATC)

BSA (10 mg/ml)

DTT (0.1 M)

dNTPs (25 mM)

PCR purification kit (Qiagen)

Platinum Taq polymerase (Invitrogen)

DNase I

2. Methods

2.1 Chromatin-IP

See protocol chromatin-IP

2.2 Amplification of CHIPed DNA

1. Assemble the following components in preparation for random primed linear amplification (Amplify IgG, AcH3, and α TF:

Per Reaction

Purified DNA 10 μ l

5X sequenase buffer (USB#70775) 4 μ l

primer A (200 μ M) 4 μ l

GTTTCCCAGTCACGATCNNNNNNNNN, HPLC purified

Total Volume 18 μ l

2. Cycling Conditions:

- a. 94 °C for 4 minutes
- b. snap cool on ice
- c. 10°C for 5 min.

3. During incubation for each reaction assemble the following components:

Per Reaction

BSA (10 mg/ml)	0.2 µl
DTT (0.1 M)	1.0 µl
dNTPs (25 mM)	0.5 µl
Sequenase (1.3 U/µl dilute stock 1:10)	<u>1.0 µl</u>
	2.7 µl

4. Add 2.6 µl of the cocktail to the reaction mixture and use the following cycling conditions:

- a. Ramped to 37° C over 9 min (1°C/20 sec) to anneal primer A
- b. 37° C for 8 min
- c. 95° C for 4 min
- d. snap cool on ice
- e. 10° C hold
- f. add 1.0 µl Sequenase (1.3 U/µl) to each sample
- g. 10° C for 5 min
- h. Ramped to 37° C over 9 min (1°C/20 sec)
- i. 37° C for 8 min
- j. repeat c through i 2X
- k. hold at 4° C

5. Purify DNA using PCR purification kit (Qiagen).

7. Amplify DNA by PCR using primer B (GTTTCCCAGTCACGATC) and Platinum Taq polymerase by assembling the following components:

Purified DNA	20 µl
10X PCR Buffer	10 µl

MgCl ₂ (25 mM)	3 μl
dNTPs (10 mM)	4 μl
Primer B	4 μl
Water	57 μl
Platinum Taq	<u>2 μl</u>
Total Reaction Volume	100 μl

- a. 95° C for 30 sec Denaturation
 - b. 40° C for 30 sec Annealing
 - c. 50° for 30 sec Annealing
 - d. 72° C for 60 sec Elongation
 - e. cycle back 29X
8. Amplified DNA is purified using Qiagen PCR purification kit.
 9. DNA size is determined by loading 0.5 μl onto 1.2% agarose gel with 100 bp ladder (expected size range 200-1000 bp).
 10. DNA was quantified using nanodrop (50-100 ng/μl or 6-9 μg total expected recovery).
 12. DNA fragmentation accomplished by DNase I treatment.

2.3. Labeling of the probe

Labeling was used with TdT and biotinylated ddATP.

2.4. Hybridization

The hybridization of the labeled DNA was performed as described in the Affymetrix® Chromatin Immunoprecipitation Assay Protocol.
(http://www.affymetrix.com/support/downloads/manuals/chromatin_immun_ChIP.pdf)